



- General Surgery :
- Dental :
- Oral Surgery :
- Otolaryngology (ENT) :
- Urology :
- Gynecology :
- Dermatology :
- Cosmetology :
- Ano Rectal Surgery :
- Ophthalmology :
- Neurosurgery :
- Gastroenterology :
- Orthopedic :
- Pulmonary/Thoracic :
- Podiatry :
- Vascular treatment :

Vitta Maxim DS15

Lumics Diode , Germany / FDA certified / C+ software



Laser Varicose Veins Clinic
Laser Piles, Fistula & Fissure
Gold Standard Treatment



Varicose Veins The different grades of piles Different types of Fistula Chronic Anal Fissure

Laser Type	Gallium-Aluminum-Arsenide Diode, Semiconductor		
Wavelength	980nm	1470nm	980nm + 1470nm
Maximum Power	15W	15WW	15W + 10W
Operation Mode	CW, Single or Repeat Pulse		
Pulse Duration	10us -3s		
Repetition Rate	0.2Hz – 50KHz		
Pilot Beam	Red Diode Laser of 650nm, Power <5mW		
Energy control	Foot switch		
Control Mode	True Color Touch Screen (7 inches, resolution 1920*600px)		
Transmission System	Medical Fibers with SMA905 Connector		
Cooling system	Internal, air and thermoelectric cooling		
Controller	Microprocessor		
output	Fiber -optic coupling		
Dimensions	23cm(W) * 18cm(L) * 13cm(H)		
Optic fiber	40um/600um/800um, 3m/400um is the standard type 400um/600um/800um, 3m/400um is the standard type		
Laser Caser	Heat resistance and water proof case		
Weight	3.0KG		

**UG-302, Dreams The Mall, Gate No. 3, L.B.S. Marg,
Bhandup (W), Mumbai - 400 078.**

Contact : 7040402020

E-mail : sales@gmerk.com Website :www.gmerk.com